



TAKE THE PRESSURE OF DOING PRESSABLE CERAMICS

CADD Anatomical Wax-Up for a press to metal restoration

CADD will design your pressable case using CADD Wax as materials. Our outsource milling center will produce the wax crown or bridge. Upon receipt of the wax restoration simply wax in the margin, invest and press.

Produce more units per day, increase profitability and take the pressure off doing pressable ceramics with CADD Anatomical Wax-up for a press to metal restoration system.

CADD Wax-up for casting in metal

Streamline your waxing department with CADD Laboratories digitally mastered rapid prototyping wax/resin copings, bridges and anatomical crowns. Minimal metal finishing required due to the ultra-smooth CADD Wax surface. Patterns are easily adjusted and burnout under normal burnout cycles. CADD Wax patterns are crafted to your design specifications then shipped to your laboratory ready to sprue, invest and cast!

CADD Wax brings a new level of consistency to your C&B department.

CALL TODAY
(661) 402-2804

